Medtronic

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9. Villavicencio RT, Rex DK, Rahmani E. Efficacy and complications of argon plasma coagulation for hematochezia related to radiation proctopathy. *Gastrointest Endosc.* 2002;55:70-74.

10. Rotondano G, Bianco MA, Marmo R, Piscopo R, Cipolletta L. Long-term outcome of argon plasma coagulation therapy for bleeding caused by chronic radiation proctopathy. *Dig Liver Dis.* 2003;35:806-810.

11. Shahbaz O, Tadros M, Chaletsky DM. Radiofrequency ablation for the management of radiation proctitis: an effective and safe treatment modality. *Gastrointestinal Endoscopy* 2014; 79(5): 566-567.

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Hemostasis of GAVE
and RP with the
HALO⁹⁰ Catheters

F. Scott Corbett, MD
Sarasota Memorial Hospital

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Comparison of tissue ablation areas using Barrx™ Ultra Long RFA Focal Catheter, Barrx™ 90 RFA Focal Catheter, Barrx™ Channel RFA Endoscopic Catheter and APC

5. Dray X, Repici A, Pedro Gonzalez P, et al. Radiofrequency ablation for the treatment of gastric antral vascular ectasia. *Endoscopy* 2014;08 [Epub ahead of print].

11. Shahbaz O, Tadros M, Chaletsky DM. Radiofrequency ablation for the management of radiation proctitis: an effective and safe treatment modality. *Gastrointestinal Endoscopy* 2014; 79(5): 566-567.

13. McGorisk T, Krishnan K, Keefer L, Komanduri S. Radiofrequency ablation for refractory gastric antral vascular ectasia (with video). *Gastrointest Endosc* 2013; 78(4): 584-8.

16. Dray X, Battaglia G, Wengrower D, et al. Radiofrequency ablation for the treatment of radiation proctitis. *Endoscopy* 2014; 46:970-6.

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▪ **Title:** RFA for the Management of Radiation Proctitis: an Effective and Safe Treatment Modality¹¹

▪ **Authors:** Shahbaz O, Tadros M, Chaletsky DM

▪ **Journal:** Gastrointestinal Endoscopy, 2014

▪ **Number of patients:** 10

▪ **Patient Characteristics:**

- All had rectal bleeding secondary to radiation proctopathy
- 6 had failed other therapies, 2 were transfusion dependent



▪ **Methods:**

- At least 1 RFA session
- Primary outcomes: resolution of bleeding, stabilization of hemoglobin

▪ **Results:**

- 9 had resolution of bleeding and the majority had improved hemoglobin
- 7 required 1 RFA session, 2 required 2 RFA sessions
- 1 failed after 5 sessions

▪ **Adverse events:** None

Images courtesy of F. Scott Corbett, MD, Sarasota Memorial Hospital, Florida Digestive Health Specialists, FL

11. Shahbaz O, Tadros M, Chaletsky DM. Radiofrequency ablation for the management of radiation proctitis: an effective and safe treatment modality. *Gastrointestinal Endoscopy* 2014; 79(5): 566-567.

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- **Title:** Treatment of chronic radiation proctopathy with RFA (with video)⁷
- **Authors:** Rustagi T, Corbett FS, Mashimo H
- **Journal:** Gastrointestinal Endoscopy, August 2014
- **Number of patients:** 39
- **Patient Characteristics:**
 - History of pelvic radiation treatment > 6 months prior
 - Recurrent hematochezia for > 3 months with endoscopic evidence of RP
- **Type of device used:** Barrx™ 90
- **Treatment method:**
 - Generally 2 pulses/site
 - Ablations placed approximately 1 mm proximal to dentate line
 - Restricted treatments to lengths ≤ 6 cm/session, and <270°/session
- **Outcomes:**
 - Mean # RFA sessions: 1.5
 - Mean follow-up: 28 months post-treatment (range 7-53 mo)
 - Mean time between sessions: 18.1 weeks
 - Complete cessation of visible bleeding in all patients
 - Overall 14% improvement in mean Hgb; 11.8 gm % → 13.5 gm %
 - Discontinuation of blood transfusion in 11/12 (92%)

ORIGINAL ARTICLE

Treatment of chronic radiation proctopathy with radiofrequency ablation (with video)

Tarun Rustagi, MD,¹ F. Scott Corbett, MD,² Hiroshi Mashimo, MD, PhD³
New Haven, Connecticut; Sarasota, Florida; Boston, Massachusetts, USA

Background: Chronic radiation proctopathy (CRP) is a common sequela occurring even many years after pelvic radiation. Current ablative therapies for bleeding ectatic vessels have the potential for deep tissue injury leading to ulcerations, perforation, and fistulas. Radiofrequency ablation (RFA) therapy avoids deep tissue injury and is a promising treatment for CRP.

Objective: To assess the long-term safety and efficacy of RFA for the treatment of CRP.

Design: Multicenter, retrospective analysis of a prospectively collected database.

Setting: Veterans Affairs Boston Healthcare System, Massachusetts; Sarasota Memorial Hospital and Suncoast Endoscopy of Sarasota, Florida.

Patients: A total of 39 consecutive patients with CRP.

Interventions: Endoscopic RFA of CRP.

Main Outcome Measurements: The primary endpoint of the study was complete resolution of rectal bleeding. Secondary endpoints included visually scored improvement of CRP on endoscopic follow-up by using a rectal telangiectasia density (RTD) grading score, improvement in hemoglobin level, and adverse events related to the procedure.

Results: A total of 39 male patients (mean [± standard deviation {SD}] age 72.9 ± 6.6 years) were included in the study. The mean number of RFA sessions was 1.49 (median 1, interquartile range [IQR] 1-2, range 1-6), with a mean interval of 18 weeks between sessions. Rectal bleeding stopped completely in all patients during the mean follow-up of 28 months (range 7-53 months). A significant improvement occurred in the mean (± SD) hemoglobin level from 11.8 ± 2 to 13.5 ± 1.6 g % ($P < .0001$). Endoscopic severity also improved significantly with an improvement in the median RTD score from 3 (IQR 2-3) to 0 (IQR 0-1) ($P < .0001$). Treatment with RFA led to discontinuation of blood transfusion and iron therapy in 92% and 82% patients, respectively.

Limitations: Retrospective analysis, lack of control group.

Conclusion: RFA therapy led to complete resolution of rectal bleeding in all treated CRP patients, with improvement in clinical and endoscopic indices without any major adverse events. Further controlled studies are needed to establish RFA as the endoscopic therapy of choice for treatment of CRP. (*Gastrointest Endosc* 2014;■:1-9.)

Adverse events:

- Transient mild to moderate anorectal pain (12%) and transient fecal incontinence (2 patients)
- One patient had a significant but hemodynamically stable anorectal hemorrhage after deviation from suggested protocol (2 ablations of 12J/cm² per treatment site), treated with hemostatic clip, no need for transfusion

7. Rustagi T, Corbett FS, Mashimo H. Treatment of chronic radiation proctopathy with radiofrequency ablation (with video). *Gastrointestinal Endoscopy* 2014; pii: S0016-5107(14)01406-0. doi: 10.1016/j.gie.2014.04.038. [Epub ahead of print.]

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
- Contraindications:
 - Pregnancy
 - Presence of gastric or colorectal ulcers
 - History of anal incontinence
 - Presence of anorectal fistulae
 - Pelvic irradiation within the last six months
- Warning:
 - Failure to discontinue platelet inhibiting and anti-thrombotic agents 7 days before and after therapy may lead to higher rates of bleeding, and the need for transfusion and/or additional intervention for hemostasis
- Most studies and case reports used the Barrx™ 90 (GAVE, RP) or Barrx™ 90 ULTRA (GAVE)
- Potential benefits of Barrx™ Channel RFA Endoscopic Catheter:
 - No repositioning of the catheter to reach target areas
 - Ease of intubation given through-the-scope catheter
 - Allows ablation within gastric folds
 - When treating RP, no need to retroflex endoscope and catheter as entire ablation zone is well visualized

12. Instructions For Use. Specific to the use of the Barrx™ Focal Catheters for the sub-indication of bleeding and non-bleeding sites in the gastrointestinal tract. Covidien, 15 Hampshire St., Mansfield, MA 02048. Part No.1062395 717-0051-01 and 717-0056-01.

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


970 Innovations and brief communications



Radiofrequency ablation for the treatment of radiation proctitis

Authors

Xavier Dray¹, Gorgio Battaglia², Dov Wengrower³, Pedro Gonzalez⁴, Alessandra Carlino⁵, Marine Camus¹, Tomer Adar³, Francisco Pérez-Roldán⁴, Philippe Marteau¹, Alessandro Repici⁵

Institutions

Institutions are listed at the end of article.

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